

Stand For Flow Cup

It is important for any flow cup to be level for the flow cup with its sample during testing of sample viscosity. Biuged instruments offers a simple type stand. It is made up of three adjustable aluminum pole, and glass plate with a spirit level .

Ordering information:

[BGD 130 -----Square Stand](#)

[BGD 131-----Tripod Stand](#)

[BGD 132---- Stainless Steel Stand](#)

[BGD 133---- Tripod Stand with double-walled thermo jacket](#)



BGD 130



BGD 131



BGD 132



BGD 133